



21  
22  
23  
24  
25  
26  
27  
28  
29  
30  
31  
32  
33  
34  
35  
36  
37  
38  
39  
40  
41  
42  
43  
44  
45  
46  
47  
48  
49  
50  
51  
52  
53  
54  
55

## Abstract

**Objective:** to explain the meaning and to illustrate technical artifacts (aliasing as well as the ring artifact) and beam hardening (metal artifact) that can be present in the dentomaxillofacial cone beam computed tomography (CBCT), and to check the accessibility of free illustrations of these artifacts in medical publications.

**Material and methods:** One observer applied five search equations using database PubMed. The exclusion criteria were: experimental studies, animal studies, studies not related to dentomaxillofacial area, and articles with closed access. There was no time limit for the search of articles. We limited our search to English and French language.

**Results:** Only 3 articles out of 434 publications were retained after application of inclusion/exclusion criteria. In these articles only 4 annotated figures were freely accessible in medical publications from PubMed.

In this paper we presented examples of aliasing, ring artifact, and beam artifacts from I-CAT, Carestream 9000 3D (Kodak), and Planmeca Promax 3D Mid CBCT. The intensity of beam hardening artifact varies from major degradation of image (i.e., subperiosteal implants, bridges, crowns, dental implants, and orthodontic fix appliances), through mean degradation (screws securing titanium mesh, head of mini-implant) to no beam hardening on metallic devices (orthodontic anchorage, orthodontic contention wire) or on dense objects (endodontic treatments, impression materials, Lego box). Some beam hardening artifacts arising from nasal piercing, hairs, or hearing aid device may be present on the image but they will not disturb the evaluation of the field of view.

**Conclusions:** reduction of aliasing artifact is related with the improvement of detectors quality. The presence of ring artifact means that CBCT device has lost its calibration. The field of view (FOV) needs to be reduce in order to avoid scanning regions susceptible to beam hardening (e.g., metallic restorations, dental implants). Finally, the accessibility to open knowledge on technique -related CBCT artifacts seems extremely limited when searching at PubMed database.

**Keywords:** CBCT, artifact, beam hardening artifact, aliasing artifact, ring artifact

56

57

## Introduction

58

59

60

61

62

63

64

65

66

67

68

69

70

71

72

The artifact represents an abnormal signal; and this is related to the conditions in whose the measurement was performed. Technical artifacts in cone beam computed tomography (CBCT) could be due to the CBCT device (aliasing, ring artifacts) and to the cone beam interaction with metallic structures inside the dentomaxillofacial area (beam hardening artifact). Artifacts may degrade the quality of the image and this could lead to a wrong diagnosis if they are not recognized and corrected. If artifacts are very strong then image obtained may be unusable and this could results with a need to repeat the exam, which is against the radioprotection principle (ALARA-as low as reasonably achievable). Therefore, dental practitioners should be aware of existence of artifacts related to the CBCT technology and understand their source. In this study we wanted to investigate the accessibility of dental practitioners to open knowledge (open access articles) on technically related CBCT artifacts. We also wanted to explain the meaning and to illustrate diverse presentations of these artifacts for clinicians to equip them with practical knowledge to perform CBCT exams better.

73

## Materials and methods

74

75

76

77

78

79

80

81

82

83

84

85

86

87

88

89

90

91

92

93

94

95

One observer performed the search using only PubMed database. The exclusion criteria were: experimental studies, animal studies, studies not related to dentomaxillofacial area, and articles with closed access. There was no time limit for the search of articles. We limited our search to English and French language. We performed 5 search equations from PubMed.

- First search equation was focused on accessible free full-text reviews about CBCT artifacts. The search equation was: `cbct[All Fields] AND ("artifacts"[MeSH Terms] OR "artifacts"[All Fields] OR "artefact"[All Fields]) AND (Review[ptyp] AND "loattrfree full text"[sb] AND "humans"[MeSH Terms])` and was performed on 05.12.2019. We found 7 articles. After full-text reading only 2 articles were accepted for this review [1, 2];
- The aim of the second search was to find open access articles on artifacts in dentomaxillofacial CBCT from a larger perspective. Our search equation was set as: `cbct[All Fields] AND ("artifacts"[MeSH Terms] OR "artifacts"[All Fields] OR "artefact"[All Fields]) AND ("humans"[MeSH Terms] AND (English[lang] OR French[lang]))` and was performed on 05.12.2019. We found 390 studies. However, after title and abstract lecture there were no articles retained for further review due to the selected exclusion criteria;
- The third search equation was related to beam hardening artifact in CBCT. Our search equation was set as: `(beam[All Fields] AND hardening[All`

96 Fields] AND CBCT[All Fields]) AND "loattrfree full text"[sb]. The search  
 97 was performed on 08.02.2020. We found 27 articles, and after applying  
 98 exclusion criteria and after full-text reading only one article was accepted  
 99 for this review [3];

- 100 • The fourth search equation was related to aliasing artifact in CBCT. Our  
 101 search equation was set as: (aliasing[All Fields] AND ("artifacts"[MeSH  
 102 Terms] OR "artifacts"[All Fields]) AND cbct[All Fields]) AND "loattrfree  
 103 full text"[sb]. The search was performed on 08.02.2020. We found 6  
 104 articles, and after applying exclusion criteria and after full-text reading we  
 105 found the same article as in third search [3];
- 106 • The fifth search equation was related to ring artifact in CBCT. Our search  
 107 equation was set as: (ring[All Fields] AND ("artifacts"[MeSH Terms] OR  
 108 "artifacts"[All Fields] OR "artifact"[All Fields]) AND cbct[All Fields])  
 109 AND "loattrfree full text"[sb]. The search was performed on 08.02.2020.  
 110 We found 4 articles. However, after title and abstract lecture there were no  
 111 articles retained for further review due to the selected exclusion criteria.

112 We stress the importance of 'annotated figures' (with arrows), because they  
 113 bring educative value. A complex radiological image presented in the  
 114 scientific article without any arrows is difficult to understand and interpret.


## 115 Results

### 116 Accessibility to open knowledge

117 We found only three articles [1-3] from a total of 434 articles that corresponded to  
 118 our search. In the study by Pauwels et al., [1] there were no figures on artifacts  
 119 available among a total number of 34 figures. In the study by Schulze et al., [2] we  
 120 found only one annotated clinical figure on aliasing artifact. Nagarajappa et al., [3]  
 121 review added one annotated clinical figure on aliasing artifact, one annotated clinical  
 122 figure of ring artifact, and one annotated clinical figure on beam hardening.  
 123 All selected articles presented very limited pictorial results. In this article we present  
 124 a wide range of artifacts (figures 1-45). All images presented below belongs to the  
 125 authors database.  
 126

### 128 Aliasing artifact

129 For each projection the voxels close to the source will be traversed by more  
 130 recorded "rays" than those close to the detector [2] because of the cone beam  
 131 divergence itself [2]. This causes aliasing which represents itself as line patterns  
 132 (moiré patterns), commonly diverging towards the periphery of the reconstructed  
 133 volume [2]. These patterns can decrease the quality of image. They can appear in  
 134 different types of dental CBCT (I-CAT, Carestream 9000 3D (Kodak)) (Figures 1,  
 135 2).  
 136  
 137



138  
139  
140  
141  
142  
143

**Fig. 1. I-CAT CBCT.** Axial view of mandible and cervical vertebra. Thick arrows: aliasing artifacts. Dashed arrow: streaks from beam hardening artifact (metal crown on tooth n°46); thin arrow: black space between metallic elements (dental implants on mandibular left side) from beam hardening artifact.




144  
145  
146  
147  
148  
149  
150

**Fig. 2. Carestream 9000 3D (Kodak).** Axial view of left mandible. Arrows: Aliasing artifact visible in soft tissues around the bone.

151  
152  
153  
154  
155  
156  
157

### Ring artifact

Ring artifacts are visible as concentric rings centered around the location of the axis of rotation [2]. They are most prominent when homogeneous media are imaged [2], such as soft tissues of the mouth floor. They are caused by defect or uncalibrated detector elements [2].




158  
159  
160  
161  
162  
163  
164  
165  
166  
167  
168  
169  
170

**Fig. 3. I-CAT CBCT.** Axial view of mandible. Arrows: ring artifact visible in the middle of mouth floor.

### Beam hardening (metal)

An X-ray beam is composed of individual photons with a range of energies [3]. As the beam passes through an object, it becomes “harder,” i.e., its mean energy increases, because the lower-energy photons are absorbed more rapidly than the higher-energy photons [3-5]. Highly absorbing materials such as metal (amalgam filling, dental implants) [3, 4] function as a filter positioned within the object [3]. Beam hardening artifacts are influenced by object density ([6]. If the emitted spectrum contains more relatively lower-energetic rays than that recorded on the

171 detector (i.e. the beam is hardened), a non-linear error (relatively too much energy  
172 recorded in the beam path behind highly absorbing materials) is induced in the  
173 recorded data [3]. In the 3D reconstruction, the error is back projected into the  
174 volume, resulting in darks streaks.[4]. Because the CBCT X-ray beam is  
175 heterochromatic and has lower mean kilovolt (peak) energy compared with  
176 conventional CT, this artifact is more pronounced on CBCT images [3]. These can  
177 be reduced using iterative reconstruction [3]. Beam hardening artifact results in dark  
178 bands or streaks between dense objects in the image [3] (Figures 5, 8, 9, 10, 11, 14,  
179 15). This artifact is related to the presence of dental amalgam (Figures 4, 5), to  
180 crowns and bridges (Figure 6-8), and to dental implants (Figure 9).  
181



182

183 **Fig. 4. I-CAT CBCT.** 3D reconstruction of mandible and maxilla. Arrows: 3D  
184 reconstruction of streaks artifacts because of beam hardening from dental  
185 amalgam fillings.

186

187

188

189


190

191

192


193

194



195  
196  
197  
198  
199  
200


**Fig. 5. Carestream 9000 3D (Kodak).** Axial view of upper maxilla. Fusion of three fields of view of axial view. Thick arrows: complex pattern of streaks from beam hardening artifact due to metallic fillings on teeth n°17, 23, 24, and 26 (dental amalgam). Dashed arrows: dark stripes between dense objects in the image (beam hardening artifact).



201  
202  
203  
204  
205  
206

**Fig. 6. I-CAT CBCT.** Axial view of mandible and of cervical spine. Thick arrows: complex pattern of streaks due to crowns and bridge present on all upper teeth. Thin arrows: aliasing artifact.






207  
208  
209  
210  
211

**Fig. 7. I-CAT CBCT.** Frontal view of mandible and of intermediary element of a bridge. Arrow: streak artefacting only a horizontal space around the metallic element. There exists no vertical streak due to beam hardening.



212  
213  
214  
215  
216


**Fig. 8. Carestream 9000 3D (Kodak).** Maxilla. A. Axial view. B. Frontal view. C. Sagittal view. Arrows showing hooks of removable partial prosthesis. Presence of streaks and black stripes due to beam hardening from metallic element.



217  
218  
219  
220  
221

**Fig. 9. I-CAT CBCT.** Axial view of maxilla. Thin arrows: Beam hardening around implants. Dashed arrows: dark stripes due to beam hardening between dense objects.

222 Especially, superiosteal implants [7] show massive beam hardening artifact  
223 making impossible the evaluation of adjacent soft and hard tissues (Figures 10-13).  
224




225 **Fig. 10. Planmeca Promax 3D Mid.** Frontal view of maxilla and mandible.  
226 Dashed arrows: upper limit of metallic artifacts from right side.  
227 Superimposition of artifact on right vestibular soft tissue, on right maxillary  
228 sinus, and on nasal fossae. Thin dashed arrows: upper limit of metallic  
229 artifact from left side. Superimposition of artifact on left vestibular soft tissue.  
230 Metallic artifacts are related to subperiosteal implants. Thin arrows: lower  
231 limit of metallic artifact that are related to the bridge. Thick arrow: metallic  
232 artifact related to the right mandibular bridge. Important superimposition of  
233 artifact on right vestibular soft tissue. Thick dashed arrow: metallic artifact  
234 related to the mandibular bridge. Few superimposition of metallic artifact on  
235 lingual side.  
236  
237  
238  
239  
240  
241  
242  
243




244  
245  
246  
247  
248  
249  
250  
251

**Fig. 11. Planmeca Promax 3D Mid.** Sagittal view of maxilla and mandible. Thin arrow: superimposition of metallic artifact from subperiosteal implants on the floor of maxillary sinus. Thick arrows: metallic artifact from subperiosteal implants superimposed on nasopharyngeal space. Thin dashed arrows: metallic artifact from subperiosteal implants superimposed on cervical spine. Thin punctiform arrow: beam hardening from mandibular dental implant.




252  
253  
254  
255  
256  
257  
258  
259

**Fig. 12. Planmeca Promax 3D Mid.** Anterior view of 3D reconstruction of maxilla and mandible. Thin arrows: metallic artifacts related to subperiosteal implants. Thick dashed arrow: massive metallic artifact related to right mandibular bridge.




260  
261  
262 **Fig. 13. Planmeca Promax 3D Mid.** 3D reconstruction of subperiosteal  
263 implant. A. Right lateral view. Thick arrow: right vestibular anterior pillar. Thin  
264 arrow: vestibular horizontal bar. Thin dashed arrow: right vestibular posterior  
265 pillar. B. Posterior view. Arrows: subperiosteal implant on palatal side. C.  
266 Left lateral view. Thin arrows: left vestibular horizontal bars. Thin dashed  
267 arrow: reinforcement bar on left palatal side.

268 Beam hardening artifact is also related to the orthodontic fix appliances (Figure 14),  
269 palatal expansion devices (Figure 15), and to orthodontic crowns (Figure 16).  
270




271  
272 **Fig. 14. I-CAT CBCT.** Axial view of mandible and orthodontic treatment,  
273 vestibular brackets, and wire. Thick dashed arrow: complex pattern of  
274 streaks (metallic artifacts) in vestibular area. Thick arrow: complex pattern of  
275 streaks (metallic artifacts) on lingual side. Thin dashed arrows: dark spaces  
276 due to beam hardening of adjacent brackets. Thin arrows: aliasing artifact.  
277  
278  
279




280  
281  
282  
283

**Fig. 15. I-CAT CBCT.** Axial view of maxillary bone. Thin dashed arrows: massive metallic artifacts around a palatal expansion device. Thin arrows: aliasing artifacts.




284  
285  
286  
287  
288  
289  
290

**Fig. 16. I-CAT CBCT.** Axial view of the mandible. Dashed arrows: orthodontic crowns on teeth n°36, and n°46. Thick arrow: Additive effect of streak artifact from two orthodontic crowns. Thin arrows: aliasing artifact.




291  
292  
293  
294

**Fig. 17. Planmeca Promax 3D Mid.** Anterior view of 3D reconstruction of maxilla and mandible. Black arrows: intermaxillary fixation screws presenting important beam hardening artifact.



295  
296  
297  
298  
299  
300

**Fig. 18. Planmeca Promax 3D Mid.** Axial view of the maxilla. Arrows: important beam hardening artifact around the head of intermaxillary fixation screw. Dashed arrow: black stripe along the long axis of the screw.



301  
302  
303  
304  
305  
306

**Fig. 19. Planmeca Promax 3D Mid.** Axial view of the maxilla. Arrows: important beam hardening artifact around the head of intermaxillary fixation screw with creation of black spaces in vestibular soft tissues surrounding the screw. Presence of black stripes along the long axis of screws. Impossibility to evaluate the area between screw and dental roots.



307  
308  
309  
310  
311  
312  
313


**Fig. 20. Planmeca Promax 3D Mid.** Frontal view of maxilla and mandible. Arrows: stripes in soft tissue around intermaxillary fixation screws due to beam hardening artifact. Double arrows: cumulation of black stripes from screws positioned at the same level on the right and the left side.

314 Some metallic elements present in the field of view may show a beam hardening  
315 artifact, however this artifact will not disturb the evaluation of dentomaxillofacial  
316 structures. Nasal piercings (Figure 21), hearing aid devices (Figures 22, 23), and  
317 metallic clips on dreadlocks (Figures 24-27) give beam hardening artifacts with  
318 streaks present outside of the dentomaxillofacial area.  
319




320  
321 **Fig. 21. I-CAT CBCT.** Axial view of maxilla. A. Arrow: piercing in left nostril  
322 wing with metallic beam hardening artifact. B. Arrow: ring of piercing in left  
323 nostril wing with metallic beam hardening artifact.  
324

325 Hearing aid device (Figures 22, 23) shows beam hardening streaks and dark space  
326 between the device and temporal bone, and streaks are present around the device and  
327 are superimposed on the skull image.  
328



329  
330 **Fig. 22. Planmeca Promax 3D Mid.** Axial view of external ear and right side  
331 of the skull. Thick arrow: hearing aid device with massive beam hardening  
332 artifact. Dashed arrow: black hole artifact close to the hearing aid device.  
333 Arrows: beam hardening artifact causing streaks superimposed on the  
334 image of the skull.






335  
336  
337  
338  
339

**Fig. 23. Planmeca Promax 3D Mid.** Sagittal view of right external ear area and of the right side of the skull. Thick arrow: hearing aid device. Dashed arrow: black hole artifact close to the hearing aid device. Arrows: beam hardening artifact causing stripes superimposed on the image of the skull.



340  
341  
342  
343  
344  
345  
346  
347  
348  
349  
350  
351

**Fig. 24. Planmeca Promax 3D Mid.** Frontal view of 3D reconstruction of mandible, maxilla, soft tissues, and hair (blue). Arrows: metallic clips on hair dreadlocks.




352  
353  
354

**Fig. 25. Planmeca Promax 3D Mid.** Frontal view of right neck area. Thick arrow: metallic clip around hair dreadlock. Thin arrow: hair dreadlock.



355  
356  
357  
358


**Fig. 26. Planmeca Promax 3D Mid.** Axial view of right neck area. Arrows: beam hardening artifact from metallic clip on the dreadlock.



359  
360  
361  
362  
363  
364  
365  
366  
367  
368  
369  
370  
371  
372  
373  
374  
375


**Fig. 27. Planmeca Promax 3D Mid.** A. Frontal view of a dreadlock. Thick arrow: metallic clip. Thin arrows: beam hardening artifact around the metallic clip. Dashed arrow: hair dreadlock. B. The 3D reconstruction of the metallic clip.

Osteosynthesis titanium orbital plates present few beam hardening artifacts giving streaks around the plate (Figures 28-30). Osteosynthesis titanium mesh provides very few beam hardening artifacts that are visible only on 2D slices (Figures 32, 33) and not present on 3D reconstruction (Figure 31). However, screws securing titanium orbital plates (Figure 29) and mesh (Figure 33) show beam hardening artifacts with presence of streaks arising in all directions around the head of the screw (Figure 29) and along the long axis of the screw (Figure 33). Osteosynthesis plates show few beam hardening artifacts that are visible only on 2D slices (Figure 35) and not present on 3D reconstruction (Figure 34).




376  
377  
378  
379  
380  
381  
382

**Fig. 28. Planmeca Promax 3D Mid.** Frontal view of the maxilla and of the orbits. Thin arrow: prebended titanium mesh for reconstruction of the floor and of the medial wall of the right orbit. Metallic artifacts from the titanium mesh and present on the right side of the mesh. Thin dashed arrow: fractured orbital bone floor fallen inside the right maxillary sinus; \*massive thickening of the right maxillary sinus mucosa.




383  
384  
385  
386  
387  
388  
389

**Fig. 29. Planmeca Promax 3D Mid.** Axial view of orbital floor. Arrows: metallic artifacts from screws positioned on the right lower orbital rim to hold the titanium mesh in place.




390  
391  
392  
393  
394  
395

**Fig. 30. Planmeca Promax 3D mid.** Anterior and right lateral view of the 3D reconstruction of the right orbit. White arrow: prebended titanium mesh. Black thin arrow: subcutaneous dense foreign body without metallic artifact (plastic?). Black thick arrow: complex fracture at the right frontozygomatic junction.




396  
397  
398  
399  
400  
401  
402

**Fig. 31. Planmeca Promax 3D mid.** Frontal view of the 3D reconstruction of frontal bone and of orbital roof. Titanium mesh fixing frontal and nasal fractures. No presence of metallic artifact.




403  
404  
405

**Fig. 32. Planmeca Promax 3D mid.** Sagittal view of frontal sinus. Arrow: titanium mesh without metallic artifact.




406  
407  
408  
409  
410  
411

**Fig. 33. Planmeca Promax 3D mid.** Axial view of frontal sinus. Short arrow: mild metallic artifact from titanium mesh superimposed on surrounding soft and bone tissue. Dashed arrows: metallic artifacts stripes from screws fixing the titanium mesh on the frontal bone. Long arrow: fracture of the right posterior wall of the frontal sinus.



412  
413  
414  
415  
416


**Fig. 34. Planmeca Promax 3D mid.** Anterior view of the 3D reconstruction of the maxilla. Black arrows: titanium osteosynthesis plates without metallic artifacts.



417  
418  
419  
420  
421  
422  
423  
424

**Fig. 35. Planmeca Promax 3D mid.** Axial view of the maxilla. Dashed arrows: metallic artifacts stripes related to beam hardening from screws fixing osteosynthesis plates to the maxilla.


Orthodontic mini-implant shows beam hardening artifact around the head of mini-implant and minor metallic artifact along the long axis of mini-implant (Figure 36).



425  
426  
427  
428  
429  
430  
431  
432  
433  
434  
435  
436  
437  
438  
439


**Fig. 36. I-CAT CBCT.** Orthodontic mini-implant. A. Axial view of the maxilla; dashed arrows: metallic artifacts stripes from, and around the head of the implant. B. Frontal view of the maxilla; arrow: mild metallic artifacts stripes from the head of the implant. C. Sagittal view of the maxilla; arrow: mild metallic artifacts stripes from the head of the implant.

Some currently used dental elements such as composite fillings (Figure 37), endodontic filling with gutta-percha (Figures 38-40), prosthetic trays (Figure 41), orthodontic contention wire (Figure 42), pediatric crown (Figure 43), and orthodontic anchorage (Bollard type) and associated screws (Figures 44, 45) do not show beam hardening artifact.




440  
441  
442

**Fig. 37. I-CAT CBCT.** Axial view of the maxilla. Arrow: composite filling on tooth n°16, without beam hardening artifact.



443  
444  
445  
446

**Fig. 38. I-CAT CBCT.** Axial view of the mandible. Roots canals sealing with gutta-percha in roots canals of teeth n°45, 44, 36 (2 canals in mesial root), and in tooth n°37. No presence of beam hardening artifact.



447  
448  
449  
450  
451  
452  
453


**Fig. 39. Carestream 9000 3D (Kodak).** Axial view of the maxilla; fusion of three fields of view. Roots canals sealing with gutta-percha in teeth n°14 (2 roots and 2 canals), 12, 23, 25, 26 (3 roots and 3 canals).






454  
455  
456  
457  
458

**Fig. 40. Carestream 9000 3D (Kodak).** Axial view of the maxilla; fusion of three fields of view. Thin arrow: root canal sealing with gutta-percha in tooth n°23. Dashed arrow: intraosseous endodontic material. Thick arrow: migration of endodontic material in subperiosteal space.




459  
460  
461  
462  
463  
464  
465  
466  
467

**Fig. 41. Carestream 9000 3D (Kodak).** A. axial view of the maxilla from a fusion of three fields of view; Thick arrows: plastic tray. Small thin arrow: dental impression material. B. Frontal view; small arrows: dental impression material. Dashed arrow: air bubble inside the dental impression material. C. Sagittal view; thick arrow: plastic tray; small arrow: dental impression material. L: lego box.




468  
469  
470  
471  
472

**Fig. 42. I-CAT CBCT.** Mandible and maxilla. A. Axial view; arrow showing upper orthodontic contention wire. B. Frontal view; arrow showing lower orthodontic contention wire. C. Sagittal view; arrows showing upper and lower orthodontic contention wires without any beam hardening artifact.




473  
474  
475  
476  
477  
478

**Fig. 43. I-CAT CBCT.** Pediatric crown on the mandible. A. Axial view; arrow: minor streak artifact around pediatric crown on deciduous tooth. B. Frontal view; arrow: minor streak artifact around pediatric crown on deciduous tooth. C. Sagittal view; arrow: no streak artifact around pediatric crown on deciduous tooth.



479  
480  
481  
482  
483  
484  
485  
486

**Fig. 44. I-CAT CBCT.** Orthodontic anchorage plate (Bollard type). A. Axial view, arrow: mild metallic artifact on the head of a screw fixing the anchorage plate. B. Frontal view; arrow: no metallic artifact on the anchorage plate. Dashed arrow: no metallic artifact on screws fixing the anchorage plate to the maxillary bone. Lower screw close to the vestibular molar root. C. Sagittal view; arrow: no metallic artifact on the anchorage plate.



487  
488  
489  
490  
491

**Fig. 45. I-CAT CBCT.** 3D reconstruction of maxilla and mandible, anterior view. Black arrows: Orthodontic anchorage plate (Bollard type) without metallic artifact.

492

## Discussion

493  
494  
495  
496  
497

The reduction of aliasing artifact is related with the improvement of detectors quality [8]. When ring artifact is found the dentist should ask for help the CBCT manufacturing company (importance of after selling service contract) as lost of calibration of CBCT can cause errors in on-screen measurements and further errors in diagnosis and treatment planning.

498  
499  
500  
501  
502  
503  
504  
505  
506

For beam hardening artifact, it is advisable to reduce the field of view (FOV) [9] to avoid scanning regions susceptible to beam hardening (e.g., metallic restorations, dental implants), which can be achieved by collimation, modification of patient positioning, or separation of the dental arches [3]. Dark spaces around implants are due to beam hardening artifact (Figure 9) and should not be interpreted as peri-implantitis. Diagnosis of peri-implantitis is clinical and could be completed using a periapical conventional intraoral X-ray [10]. The zirconia implants produce more artifacts and its images are more affected by the different protocols than titanium implants [11].

507  
508  
509  
510  
511  
512  
513

The importance of beam hardening artifact is related to the density of material and may vary from major perturbation of image with subperiosteal implants (Figures 10-13) to no perturbation with metallic orthodontic anchorage (Figures 44, 45). Beam hardening artifacts may vary according to the type and to the alloy material used in manufacturing of osteosynthesis screws. Screws used to fix orbit titanium plate (Figures 28, 29) shows less artifact than screws to fix titanium mesh on frontal bone (Figures 31-33). This may be of importance when checking if screws are not in

514 contact or damage roots of teeth after trauma or orthognathic surgery. With presence  
 515 of important beam hardening artifact around screws it becomes impossible to clearly  
 516 evaluate a close relationship between screws and dental roots.  
 517 A check list procedure should be installed in dental practice to avoid unnecessary  
 518 beam hardening artifacts. Especially, glasses, necklace, removable prostheses,  
 519 piercings (whose can be removed), earrings, hair clips should be removed before  
 520 CBCT examination. Hearing aid device could remain in place as it does not artifact  
 521 the dentomaxillofacial area of interest, and allows the contact with the patient.  
 522 Moreover, beam hardening artifacts related to hearing aid devices will be visible  
 523 only in large fields of view in dentomaxillofacial radiology. Dreadlocks clips give  
 524 beam hardening artifact outside of dentomaxillofacial area, and can stay in place  
 525 during the scanning.

526  
 527 Finally, when looking at PubMed database the accessibility to open knowledge on  
 528 technically related CBCT artifacts seems extremely limited with only 4 figures  
 529 freely accessible; thus the need to share practical knowledge related to CBCT  
 530 images. The author fully supports values of open science – sharing knowledge has  
 531 the potential to increase the quality of science system. All images are in CC-BY-SA  
 532 license and they can be freely reproduced, under condition of citing correctly  
 533 the source.

- 535 • **Acknowledgements:** none
- 536 • **Funding sources statement:** this study does not received any funding
- 537 • **Competing interests:** Prof R. Olszewski is Editor-in-Chief of Nemesis
- 538 • **Ethical approval:** There was no need for ethical committee approval for this  
 539 study
- 540 • **Informed consent:** There was no need for informed consent for this study

541 **Authors contribution:**

Author	Contributor role
Olszewski Raphael	Conceptualization, Data curation, Investigation, Methodology, Resources, Validation, Writing original draft preparation, Supervision, Writing review and editing

542 **References**

543 1. Pauwels R, Araki K, Siewerdsen JH, Thongvigitmanee SS. Technical aspects of  
 544 dental CBCT: state of the art. Dentomaxillofac Radiol 2015;44:20140224.  
 545

- 546 2. Schulze R, Heil U, Gross D, Bruellmann DD, Dranischnikow E, Schwanecke U,  
547 Schoemer E. Artefacts in CBCT: a review. *Dentomaxillofac Radiol* 2011;40:265-  
548 273.
- 549  
550 3. Nagarajappa AK, Dwivedi N, Tiwari R. Artifacts: The downturn of CBCT image.  
551 *J Int Soc Prev Community Dent* 2015;5:440-445.
- 552  
553 4. Esmaeili F, Johari M, Haddadi P, Vatankhah M. Beam hardening artifacts: Com-  
554 parison between two cone beam computed tomography scanners. *J Dent Res Dent*  
555 *Clin Dent Prospects* 2012;6:49-53.
- 556  
557 5. Ibraheem I. Reduction of artifacts in dental cone beam CT images to improve the  
558 three dimensional image reconstruction. *J Biomed Sci Eng* 2012;5:409-415.
- 559  
560 6. Nabha W, Hong YM, Cho JH, Hwang HS. Assessment of metal artifacts in three-  
561 dimensional dental surface models derived by cone-beam computed tomography.  
562 *Korean J Orthod* 2014;44:229-235.
- 563  
564 7. Streel R. Oral rehabilitation with implant dentures. *Dent News* 1987;9:36, 39-41  
565 contd.
- 566  
567 8. Xie H, Tang X. Optimization of data acquisition in axial CT under the framework  
568 of sampling on lattice for suppression of aliasing artifacts with algorithmic detector  
569 interlacing. *Med Phys* 2017;44:6239-6250.
- 570  
571 9. Shokri A, Jamalpour MR, Khavid A, Mohseni Z, Sadeghi M. Effect of exposure  
572 parameters of cone beam computed tomography on metal artifact reduction around  
573 the dental implants in various bone densities. *BMC Med Imaging* 2019;19:34.
- 574  
575 10. Steiger-Ronay V, Krcmaric Z, Schmidlin PR, Sahrman P, Wiedemeier DB,  
576 Benic GI. Assessment of peri-implant defects at titanium and zirconium dioxide im-  
577 plants by means of periapical radiographs and cone beam computed tomography: An  
578 in-vitro examination. *Clin Oral Implants Res* 2018;29:1195-1201.
- 579  
580 11. Vasconcelos TV, Leandro Nascimento EH, Bechara BB, Freitas DQ, Noujeim  
581 M. Influence of cone beam Computed Tomography settings on implant artifact pro-  
582 duction: zirconia and titanium. *Int J Oral Maxillofac Implants* 2019;34:1114-1120.
- 583